

DrugConfirm Advanced Cup CLIA-WAIVED



drugconfirm™ advanced
instant multi-drug test kit

Screen for the most crucial drugs of abuse today with DrugConfirm Advanced. This self-contained cup lets you collect samples easily, read lab-accurate results in minutes, and minimize exposure to urine. Test for multiple drugs in one sample, so you can take out the guesswork and feel confident in the results.

PRODUCT FEATURES:

- Can detect up to 14 common drugs of abuse plus adulterants
- Multiple configurations available
- Results within minutes
- CLIA-waived, FDA 510(K) & CE Marked
- Applicable for use across industries from pre-employment testing and government, to clinical and treatment centers
- 25 cups per box

TESTS FOR:

- Amphetamine (AMP)
- Barbiturates (BAR)
- Benzodiazepines (BZO)
- Buprenorphine (BUP)
- Cocaine (COC)
- Ecstasy/Molly (MDMA)
- Marijuana (THC)
- Methadone (MTD)
- Methadone Metabolite (EDDP)
- Methamphetamine (mAMP)
- Morphine/Opiates (MOR/OPI)
- Oxycodone (OXY)
- Phencyclidine (PCP)
- Propoxyphene (PPX)
- Tricyclic Antidepressants (TCA)



CHOOSE DRUGCONFIRM ADVANCED FOR YOUR TESTING NEEDS

The DrugConfirm Advanced Drug Test Cup is a fully integrated, self-contained drug test kit for detecting drugs and their metabolites in urine. Our multi-panel drug screen eliminates the possibility of donor tampering if an adulteration test is integrated and provides an easy solution for on-site testing that avoids the handling of urine.

INSTRUCTIONS

1. Tear the foil bag open, remove test cup. Label the device with donor information.
2. Open test cup lid. Urinate directly into the test cup. Be sure to fill the test cup with urine sample between a minimum of 30 ml to a maximum of 110 ml. The minimum and maximum levels are marked on the cup.
3. After urine specimen has been collected, close the lid securely and return cup to specimen collector.
4. For specimen collector: Put on gloves provided before receiving the cup. Peel off label to reveal test result. Read test result at 5 minutes. Results expire after 10 minutes.

INTERPRETATION OF RESULTS

NEGATIVE:

Two color lines appear. One color line in the Control (C) region and another line next to the test (T) region.




POSITIVE:

One color line appears in the control (C) region and no line appears in the test (T) region.



INVALID:

Control line fails to appear. In case of invalid results, review procedure and repeat test using a new test device. If problem persists, discontinue using the lot immediately and contact your supplier.



IMPORTANT INFORMATION

Do not use after the expiration date.

The test device should remain in the sealed pouch until use.

Results are preliminary and positive results do not always mean that drugs are present in the donor. Certain foods and medicines, diet plan drugs, and nutritional supplements may yield a false positive test result with this product. Similarly, if the sample is diluted or the sample is adulterated, the product may yield negative results and not detect a particular drug. Laboratory tests confirm the preliminary results.

TEST	DESCRIPTION	DESCRIPTION
AMP	AMPHETAMINE	1000 NG/ML
BAR	BARBITURATES	300 NG/ML
BZO	BENZODIAZEPINE	300 NG/ML
BUP	BUPRENORPHINE	10 NG/ML
COC	COCAINE	300 NG/ML
MDMA	ECSTASY	500 NG/ML
THC	MARIJUANA	50 NG/ML
MTD	METHADONE	300 NG/ML
EDDP	METHADONE METABOLITE	300 NG/ML
mAMP	METHAMPHETAMINE	1000 NG/ML
MOR/OPI	MORPHINE/OPIATES	300/2000 NG/ML
OXY	OXYCODONE	100 NG/ML
PCP	PHENCYCLIDINE	25 NG/ML
PPX	PROPOXYPHENE	300 NG/ML
TCA	TRICYCLIC ANTIDEPRESSANTS	1000 NG/ML

*For forensic use only when EtG or K2 are included