



Omixon Holotype HLA™

Make the switch to NGS for HLA

January, 2016



Peter Meintjes, PhD



Omixon Holotype HLA

- Believe the Data
 - 2000+ routine clinical samples at CHOP
 - 253 samples for ASHI validation study
 - 2500+ samples in worldwide EAP
 - 200 samples in CE-IVD Performance Evaluation Study
- Technician friendly protocol
 - All loci pooled together
 - Only 1 bead cleanup
 - 4 hours hands on time and 3.5 day TAT
 - Automation
- Scalability
 - 8 - 192 samples per run
 - 1, 5, 7 or 11 loci
- Entire team focused on NGS for HLA

Omixon Manufacture & Supply



- US-based contracted manufacturing
- ISO-13485 and GMP
- Operating Statistics
 - 83,000 active customers
 - 2,800 orders shipped per day
 - 20,000 SKUs under management
 - 100+ OEM product lines
 - 4,500+ finished goods produced per month
- Worldwide distribution



Omixon Support and Services



- Presales
 - Comprehensive NGS Cost Analysis
 - Protocol viewing at an Omixon Reference Lab
 - Hands on workshop at an Omixon Reference Lab
 - Site Visit Consult/Presentation
- Site Preparation
 - Full equipment, reagent and consumable check
- Evaluation and Training
 - 1 week onsite training
 - Data analysis walkthrough
 - Ongoing remote and onsite support
- Validation
 - 200 samples and 50 blind
 - Complex case dataset
- Clinical Production
 - Parallel testing

EAP Results



- 24 labs, 2530 samples, 26286 alleles
- 8788 alleles with known typings
- 1.7% background FTA (failure to amplify) rate

Locus	Total alleles (with known types)*	Unique alleles
HLA-A	1490	56
HLA-B	1468	84
HLA-C	1481	47
HLA-DQB1	1491	23
HLA-DRB1	1530	53

- * = some known types were missing

Concordance



Category	Count	Percent
Alleles with known typings	8878	100%
Total discordances	333	3.8%
Incorrect known typing	147	1.7%
Discordant with QC failure	85	1%
Discordance due to novelty	63	0.7%
Still under investigation	35	0.4%

- Concordance = 2 field concordance

Result Statistics



Locus	HLA-A	HLA-B	HLA-C	HLA-DQB1	HLA-DRB1
Typed alleles	1490	1468	1481	1491	1530
Sensitivity	98.65%	97.48%	98.18%	95.71%	96.41%
Specificity	99.98%	99.97%	99.97%	99.88%	99.96%
PPV	98.66%	97.57%	98.18%	95.71%	97.04%
NPV	99.98%	99.97%	99.97%	99.88%	99.95%
TCC (Accuracy)	99.96%	99.95%	99.94%	99.76%	99.91%

- PPV = Positive Predictive Value
- NPV = Negative Predictive Value
- TCC = Type Correctly Classified (TP + TN / N) = Accuracy

Questions?



Peter Meintjes, PhD

Chief Commercial Officer

peter.meintjes@omixon.com

+1 (203) 947 2772

OMIXON.COM

