


hsTnT

As a five level, liquid-stable product, VALIDATE® hsTnT is tailor-made to specifically challenge the reportable range of Roche's Troponin T Gen 5 STAT (TnTG5) assay, eliminating the need for manual dilutions.

This FDA listed material allows clinical laboratories to fulfill their Linearity and Calibration Verification requirements for high sensitivity Troponin T.



Results of Lot 45AA16918R on Roche Cobas® TNTG5 assay. Mean versus target regression: $y = 1.032x - 2.724$ ($r^2 = 0.9997$). Limits applied as 50% of the Total Allowable Error (the greater of 2 ng/L or 22% of recovered means).

VALIDATE® hsTnT

Order No. 405ro

Typical Recovery: 12 - 936 ng/L
 Matrix: Human Serum
 Form: 1.5 mL, Liquid-Frozen
 Storage: -10° to -25°C
 Open-Vial: 4 freeze-thaws through expiration

Units: ng/L	L1	L2	L3	L4	L5
Target	10.0	223.0	436.0	649.0	862.0
Mean	10.0	223.0	445.3	675.3	883.3
% Diff	0.0%	0.0%	9.3%	26.3%	21.3%

EASY • FAST • EFFICIENT



LIQUID, READY-TO-USE

Increase productivity and avoid the time, error and inconvenience of reconstitution. Supplied in dropper bottles, simply bring to room temperature and dispense into the sample cup, reducing the risk of contamination.



MULTI-USE & OPEN-VIAL STABILITY

Multiple freeze-thaws, through to the expiration of the product, provide for 2+ testing cycles. Cap the bottle, return to storage and use again.



TROUBLESHOOTING

As a multi-use product, VALIDATE® is readily accessible for use after preventative maintenance or to troubleshoot across the full reportable range.



LINEARITY & PEER GROUPS

5 Level products with free data reduction services offers instant linearity results via MSDR® with ready access to Peer Statistics.



QUALITY

Products are U.S. FDA cleared or listed and are CE marked. Our facility is registered with the U.S. FDA and our QMS is certified to ISO 13485.



ACCREDITATION

For use at installation, or during the laboratory's routine 6-month testing interval, rely on VALIDATE® to satisfy ISO 15189, CAP, JCI and local regulatory requirements, including AMR validation.

