

Midbrook Device Washer

Multiple mechanical actions to clean and disinfect multiple surgical device styles.



Cleans and disinfects:

Canulated Instruments
Endoscopy Devices

Lumen Instruments
Orthopedic Devices

Powerful cleaning and disinfection for all surgical instruments

For use in:

Central Processing

Ambulatory Surgical Centers

Operating Rooms



- Fully automated
- Small footprint
- Flexibility
- Efficient disinfection
- Ozonated cold water flush
- Enzymatic, agitated submersion
- Heated submersion
- Heated microbubble fluid flush
- Final external/internal rinse
- Filtered air blowout

Assure consistent, repeatable results that hand cleaning or ultrasonics alone cannot duplicate