



Rx Only

For Topical Dermatological Use Only

Product Description

EpiCeram® Skin Barrier Emulsion is a steroid-free, fragrance-free, ceramide-dominant formulation.

Indications for Use

EpiCeram® Skin Barrier Emulsion is to be used to treat dry skin conditions and to manage and relieve the burning and itching associated with various types of dermatoses, including atopic dermatitis, irritant contact dermatitis, and radiation dermatitis. EpiCeram® Skin Barrier Emulsion helps to relieve dry, waxy skin by maintaining a moist wound and skin environment, which is beneficial to the healing process.

Contraindications

EpiCeram® Skin Barrier Emulsion is contraindicated in persons with known hypersensitivity to any of the components of the formulation.

Warnings

EpiCeram® Skin Barrier Emulsion does not contain a sunscreen and should always be used in conjunction with a sunscreen in sun exposed areas. In radiation dermatitis and/or in conjunction with ongoing radiation therapy apply following radiation therapy. Do not apply within 4 hours prior to radiation therapy. Apply twice daily or as indicated by the radiation therapist. After application, a temporary tingling sensation may occur (10 to 15 minutes). Keep this and similar products out of the reach of children. Follow directions for use. If condition does not improve within 10 to 14 days, consult a physician.

Precautions and Observations

For the treatment of any dermal wound, consult a physician.

- Use EpiCeram® Skin Barrier Emulsion only as directed.
- EpiCeram® Skin Barrier Emulsion is non-toxic, however it is for external use only and should not be ingested or taken internally.
- If clinical signs of infection are present, appropriate treatment should be initiated. If clinically indicated, use of EpiCeram® Skin Barrier Emulsion may be continued during the anti-infective therapy.
- If the condition does not improve within 10-14 days, consult a physician.

- EpiCeram® Skin Barrier Emulsion does not contain a sunscreen and should always be used in conjunction with a sunscreen in sun exposed areas.
- In radiation dermatitis and/or in conjunction with ongoing radiation therapy, apply following radiation therapy.
- Do not apply within 4 hours prior to radiation therapy.
- Apply twice daily or as indicated by the radiation therapist.
- Following the application of EpiCeram® Skin Barrier Emulsion a temporary tingling sensation may occur (10-15 minutes).
- Keep this and other similar products out of the reach of children.

Instructions for Use

Apply in a thin layer to the affected skin areas 2 times per day (or as needed) and massage gently into the skin. If the skin is broken, cover EpiCeram® Skin Barrier Emulsion with a dressing of choice.

Ingredients

Capric Acid, Cholesterol, Citric Acid, Conjugated Linoleic Acid, Dimethicone, Disodium EDTA, E. Cerifera (Candelilla) Wax, Food Starch Modified Corn Syrup Solids, Glycerin, Glyceryl Stearate, Hydroxypropyl Bispalmitamide MEA (Ceramide), Palmitic Acid, PEG-100 Stearate, Petrolatum, Phenoxyethanol, Potassium Hydroxide, Purified Water, Sorbic Acid, Squalane, Xanthan Gum.

How Supplied

EpiCeram® Skin Barrier Emulsion are supplied as follows:
NDC 51013-800-90: 90 gram tube
NDC 51013-800-50: 50 gram tube

Store at 15°C to 30°C (59°F to 86°F). Do not freeze.

EpiCeram is a registered trademark of PuraCap Pharmaceutical LLC, South Plainfield, NJ 07080

Distributed by:
PuraCap Pharmaceutical LLC, South Plainfield, NJ 07080
www.EpiCeram-US.com

U.S. Patent No. 5,643,899.

**Rx only — Prescription Medical Device;
Federal Law restricts this device to sale
by or on the order of a physician.**



P190-00
Issued June 2011