

Puritan[®] Fecal Opti-Swab[®] Collection & Transport System.

Optimize your specimen collection and preservation of fecal enteric bacteria with Puritan[®] Fecal Opti-Swab[®]. A liquid base ready-to-use system for rectal swabbing or direct collection from fecal sample.




HOW IT WORKS

Simple to use – eliminates unnecessary specimen handling steps.

- Our HydraFlock[®] swab assures superior collection and elution of sample
- Only a small amount of sample is needed, making handling and transporting easy
- Versatile system allows for use with traditional plating methods or liquid based automation platforms

TESTING APPLICATIONS

While different types of bacteria can cause enteric infections routine stool cultures can generally be used to screen for public health issues such as *Salmonella spp.*, *Shigella spp.*, and *Campylobacter spp.* Can also be used to culture: *E.coli*, *V. parahaemoliticus*, *Escherichia coli* O157:H, *Y. enterocolitica*, *Shigella flexneri*, *C. jejuni*, *S. typhimurium*, *C. difficile*, and *E. faecalis* (VRE).



Optimize your specimen collection.



CONVENIENT, ALL-IN-ONE KIT.

Each System is enclosed in a sterile, easy-peel pouch.

LEAK-PROOF, SELF-CENTERING CAP

Self-centering cap is designed to facilitate easy swab removal and disposal.

LIQUID MEDIUM

2ml Cary-Blair medium is compatible with manual and automated specimen processors.

MEDICAL GRADE MATERIALS

The vial allows for safe collection and transport of specimens. Made in the US.

MAX FILL

Indicator line showing maximum sample quantity.



SUPERIOR ABSORPTION

Puritan® HydraFlock® provides superior absorption, outperforming competition.

MOLDED BREAKPOINT

Handle is sturdy and flexible. Offers a molded breakpoint with printed indicator line for safe insertion.



PURITAN® EXCLUSIVE

HydraFlock® is a patented proprietary manufacturing design process offered only by Puritan®.

This creates a three-dimensional microstructure fiber that enhances collection and release capabilities.

PURITAN® FECAL OPTI-SWAB® COLLECTION & TRANSPORT SYSTEM

Item No.	Product Description	Volume	Tube length x diameter	Packaging pack, weight
CB-206	Media: Cary-Blair medium Swab: One flock-tipped collection swab (HydraFlock 3506-H)	2 ml	12 x 80 mm	6/50, 9 lbs



207-876-3311 • puritanmedproducts.com
sales@puritanmedproducts.com
 Puritan Medical Products Co.
 31 School Street, Guilford, Maine 04443-0149 USA
 ISO 9001:2008 ISO 13485:2003 CE



- facebook.com/PuritanUSA
- twitter.com/PuritanUSA
- linkedin.com/company/PuritanUSA
- youtube.com/PuritanUSA