


# Introduction – The ABC's of DNA


DNA (deoxyribonucleic acid) is a negatively-charged doublestranded chain of nucleotides.

- Purine Bases
  - Adenine (A)
  - Guanine (G)



- Pyrimidine Bases Thymine (T)

  - Cytosine (C)



Copyright © 2013 Crosscutting Concepts, LLC. All Rights Reserved. www.CrosscuttingConcepts.com

#### Introduction – The ABC's of DNA

This chain also includes a 'backbone' made from molecules of sugar (deoxyribose) and phosphate.

In 1950, biochemist Erwin Chargaff discovered that the number of Adenine and Thymine were equal and that the same was true for Guanine and Cytosine.

What conclusions can you draw from Chargaff's discovery?

Copyright © 2013 Crosscutting Concepts, LLC. All Rights Reserved. www.CrosscuttingConcepts.com



# **DNA Double Helix**

DNA has directionality – each chain has a 'top' and a 'bottom' end.

- 'Top' end is called **Five Prime** (5')
- 'Bottom' end is called **Three Prime** (3')

Two nucleotide chains together wind into a **double helix.** 

Hydrogen bonds between paired bases hold the two DNA strands together..

DNA strands are **antiparallel** and complementary to each other.


5' ACGTC 3' 3' TGCAG 5'

Copyright © 2013 Crosscutting Concepts, LLC. All Rights Reserved. www.CrosscuttingConcepts.com

# **Mitochondrial DNA** (mtDNA) is different from the DNA found in the nucleus of the cell.

- Ring shaped
- Less Protection from mutation
- Codes for only 37 genes
- Has about 16,000 bp
- mtDNA is haploid not diploid


Because the tiny amount of mitochondria in human sperm is destroyed during fetal development, mtDNA is always inherited from the mother's egg – a male's mtDNA comes from his mother, but will not pass on to his offspring.



**Mitochondrial DNA** 

A lane is the straight line traveled by the DNA molecules between the well and the positive end of the gel.


A size ladder is a solution of DNA molecules of known lengths which is usually run on the first and/or last lanes of the gel to let us measure the sizes of the DNA pieces being tested.



# **PCR – Preliminary Step**

In order to perform PCR, a number of 'ingredients' and tools are necessary. The most important ones are:

- DNA sample
- **Polymerase** Polymerase binds to open DNA and synthesizes a complementary strand from primer bound to the original strand; Taq polymerase is a special polymerase which does not degrade at high temperatures
- Primers Short sequences of DNA that bind to specific complementary (matching) regions of the DNA sample – called the target sequence – that we want to amplify
- Nucleotides Raw material for the polymerase to synthesize new DNA strands
- **Thermocycler** Machine used to precisely raise and lower the temperature of the DNA samples for specific lengths of time



www.CrosscuttingConcepts.com



To access the complete lecture presentation, login to the Downloads area at <u>www.CrosscuttingConcepts.com</u>

# **An Identity Crisis<sup>™</sup> DNA Typing**

# crosscuttingconcepts.com

An Identity Cris

