




SpecuLume[®]


DISPOSABLE CHEMILUMINATED
VAGINAL SPECULUM




EASY TO USE



Remove speculum
and insert from
packaging



Bend to activate
light insert and
shake to further
activate



Slide light insert
firmly into distal end
of upper blade

CHEMILUMINESCENT INSERT PROVIDES OPTIMAL ILLUMINATION
FOR **OVER 15 MINUTES**

REQUIRES **NO ELECTRICITY OR BATTERIES**

SINGLE-USE MITIGATES RISK OF CROSS-CONTAMINATION

UNIQUE LIGHT SOURCE GENERATES **NO HEAT**

SIMPLY DISCARD THE ENTIRE PRODUCT AS CLINICAL WASTE
AFTER USE

ELIMINATES THE COST OF REPROCESSING AND CLEANING

AVAILABLE IN **THREE SIZES**: SMALL, MEDIUM, LARGE

FOR MORE INFORMATION OR TO REQUEST A SAMPLE,
VISIT CYALUMEMED.COM