



COMING 1st Quarter 2013

CONTACT


Sharon Krout
717.235.6807 x 2303
skrout@oakworks.com
www.oakworks.com

Preliminary Product Specifications

•	1 3
STANDARD	
Motion I	19" - 35" Height Range
Motion 2	I5°/5° Trendelenburg/Reverse Trendelenburg
Motion 3	80° Backrest/Fowler
Motion 4	90° Leg Rest
Foot Control	Supports all powered motions
Head Rest	10" long x 10" wide, manually adjustable
Paper Roll Holder	Accommodates 18" rolls (rolls not included)
Table Voltage Options	110V/60Hz or 230V/50Hz
Patient Capacity	600 lbs.
3 Section Top	75" long x 27" wide
Backrest/Fowler	26" long x 27" wide
Seat Section	23" long x 27" wide
Leg Rest Section	12" long x 23" wide
Padding	3" thick multi-layered
Upholstery	Premium, PVC free TerraTouch™
Enclosure	Black
Table Weight	300 lbs.
Regulatory	Pending
OPTIONS	
Arm Rests	Auto adjusts to backrest angle, swing out 90°, 24" L x 5" W, removable
Retractable Foot Rests (Stirrups)	Pull out from under seat section, fully adjustable
Procedure Tray	12.25" x 10.25", 2 Quarts/1.9 Liters
Enclosure Upgrade	Designer Burl Wood finish
ACCESSORIES	
T-rail Adapter	Attaches to the chair for any US format T-rail attachable accessory
Memory Positions	Hand control required
	•



Top Dimensions



Oakworks® 300 Series Procedure Chair Line

Model

300	Base Model
301	Base Model w/ Arm Rests
320	Base Model w/ Retractable Foot Rests (Stirrups)
321	Base Model w/ Arm Rests & Retractable Foot Rests (Stirrups)

Model

	300A	Aesthetic Upgrade Model
	301A	Aesthetic Upgrade Model w/ Arm Rests
	320A	Aesthetic Upgrade Model w/ Retractable Foot Rests (Stirrups)
	321A	Aesthetic Upgrade Model w/ Arm Rests & Retractable Foot Rests (Stirrups)

Add "FP" to above model for optional Procedure Tray.

