

3DEEP RF on TV and in Magazines!



### How can 3DEEP RF help me?

3DEEP RF is the most effective RF technology available today for full body skin tightening and body contouring. 3DEEP RF is ideal for tightening lax skin, improving cellulite appearance, and reducing body circumferences. Benefit from a safe, totally painless technology that delivers precise energy to target all treatment areas.

- Tighten lax skin
- Wrinkle reduction
- Improvement of acne scars
- Reduce circumferences and cellulite
- Non-invasive treatment of hyperhidrosis

Discover 3DEEP® RF –  
Your skin, only better!

Schedule a consultation today!

*Intensif™*



Impressive  
Results  
with minimal  
pain and  
downtime

International (non-US) patient brochure

 **EndyMed™**  
THE 3DEEP® COMPANY

[ e ] [info@endymed.com](mailto:info@endymed.com) [ w ] [www.endymed.com](http://www.endymed.com)  
© 2013, EndyMed Medical Ltd. All rights reserved.

 Like us on Facebook at [www.facebook.com/EndyMed](http://www.facebook.com/EndyMed)

StudioGE.com

PB6003-?



 **EndyMed™**  
THE 3DEEP® COMPANY

# Welcome to the next level in Radiofrequency skin treatments – EndyMed's new RF Microneedle Treatment Handpiece



Introducing the Intensif, the most technologically advanced RF treatment therapy available today. The Intensif treatment handpiece provides the most effective dermal skin remodeling, deep within the skin to revive and rejuvenate the natural production of new collagen and elastin. Treatment with the Intensif triggers the skin's natural healing process, rejuvenating and refreshing the skin while reducing the appearance of wrinkles, acne and traumatic scars in all skin types.

And the benefits are long lasting – with the Intensif, your skin will continue to improve up to a year after your last treatment session!

## The Intensif Advantage:

1. The most effective wrinkle reduction and acne scar treatment that offers consistent results with minimal pain and downtime. The perfect option for the treatment of acne scars, deep wrinkle, stretch marks and dilated pores.
2. The Intensif's sterile, high quality Gold Plated Microneedles are tiny (300microns thin) and tapered for smooth insertion to ensure minimal pain and downtime.
3. The Intensif's shock-free motion needle insertion mechanism allows for painless, effective and consistent treatments every time.
4. The high precision needle depth control offered by the Intensif provides for the most effective treatment strength suitable to the particular treatment area, for any skin type.

## 3DEEP RF Yields Visible Results

### Wrinkle Reduction / Skin Tightening



Before treatment



After 1 treatment

### Acne Scars




Before treatment



1 month after 2 treatments

## How it works

The EndyMed Intensif Microneedles emit heat into the dermis, triggering the production of new collagen



Photos 1: Courtesy of Itzhak Vider, M.D., Hertzlia, Israel

Photos 2: Courtesy of Monica Elamn, M.D., Rishon Le'Zion, Israel